

USP <644> Conductivity PT EXPRESS - Available year-round

Art. ID USP-PHPT-006B
Unit each

Description

Quantitative Proficiency Test Express for laboratories performing Conductivity by USP <644> or equivalent. Participants report quantitative results for conductivity, testing for material's conductivity, in micro siemens per centimeter (uS/cm), for their sample(s). Samples provided two vials containing 10g of material. Participants are asked to determine the conductivity of target analyte in the sample per the USP monograph or equivalent. This product is available year-round. It does not need to be ordered during an open PT Study. /// Target Analyte: Erythriol, powdered or granulated solid

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
Target Analyte, powdered or granulated solid	Erythriol					