

USP <197U> PT EXPRESS for Identification by UV/Visible Spectroscopy - Available year-round

Art. ID USP-PHPT-010B
Unit each

Description

Qualitative Proficiency Test Express for laboratories performing UV/Visible Spectroscopy by USP <197U> or equivalent. Participants report qualitative results for their samples. Samples provided as a three sample panel. One of three samples will be positive for target analyte. Participants are asked to identify the positive sample by UV/Vis ID test per USP monograph. This product is available year-round. It does not need to be ordered during an open PT Study. /// Positive Sample: Sorbic Acid, powdered or granulated solid /// Negative Samples: Tartaric Acid, powdered or granulated solid /// Sorbic Acid Reference Material Standard will be provided with the Proficiency Test Samples

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
Positive Sample, powder ed or granulated solid	Acetaminophen	[103-90-2]				
Negative Samples, powde red or granulated solid	Caffeine	[58-08-2]				