

## **USP <921> Water Determination Proficiency Test**

Art. ID USP-PHPT-011

Unit each

### **Description**

Quantitative Proficiency Test (PT) for laboratories performing Water Determination testing by USP <921> or equivalent. Participants report quantitative results for their liquid sample(s). Samples provided two vials containing 10mL of material. Participants are asked to determine the percent (%) water content in a liquid pharmaceutical product sample per the USP monograph or equivalent. /// Sample Material: Isopropyl Alcohol, between 60-99% (w/w)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
Sample Material: Isopro pyl Alcohol, between 60 -99% (w/w)	2-Propanol	[67-63-0]				