

USP <731> Loss on Drying Proficiency Test

Art. ID USP-PHPT-013

Unit each

Description

Quantitative Proficiency Test (PT) for laboratories performing Loss on Drying testing by USP <731> or equivalent. Participants report quantitative results for their sample(s). Samples provided in duplicate. Participant will receive two vials containing 2g of material. Participants are asked to determine the percent (%) loss on drying for the sample per the USP monograph or equivalent. /// Sample Material: Corn Starch, powdered solid

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
Sample Material, powder ed solid	Corn Starch					