

USP <831> Proficiency Test for Refractive Index

Art. ID USP-PHPT-014

Unit each

Description

Qualitative Proficiency Test (PT) for laboratories performing Refractive Index tests by USP <831> or equivalent. Participants report qualitative results for their samples. Samples provided as a three sample panel. One of three samples will be positive for target analyte. Participants are asked to identify the positive sample by Refractive Index ID test per USP monograph. /// Positive Sample: Isopropyl Alcohol, 1.5mL /// Negative Samples: 1-Propanol, 1.5mL /// Isopropyl Alcohol Reference Material Standard will be provided with the Proficiency Test Samples

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
Positive Sample: Isopro pyl Alcohol, 1,5mL	2-Propanol	[67-63-0]				
Negative Samples: 1-Pro panol, 1,5mL	1-Propanol	[71-23-8]				