

4-HYDROXYHIPPURIC ACID (RING-13C6, 99%), 100 µg/mL in Methanol

Art. ID CIL-CLM-11190-1.2

Unit 1,2 mL

Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	4-HYDROXYHIPPURIC AC [n/a] (RING-13C6, 99%)		100	µg/mL		