

2-METHYL-3-PHENYLBENZOIC ACID (PHENYL-13C6, 98%), 100 µg/mL in Acetonitrile

Art. ID CIL-CLM-11214-1.2
Unit 1,2 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1648 II E2

Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2-METHYL-3-PHENYLBEN [n/a] C ACID (PHENYL-13C6, 98 %)		100	µg/mL		