

(+/-)-TRANS-3'-HYDROXYCOTININE (METHYL-D3, 98%) (95% chemical purity), 100 µg/mL in Methanol

Art. ID CIL-DLM-4787-1.2
Unit 1,2 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

Description

Chemical purity: 95%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	trans-3'-Hydroxy Cotinine-d3	[159956-78-2]	100	µg/mL		