

METHYLPHOSPHONIC ACID (METHYL-D3, 98%), 1000 µg/mL in Methanol

Art. ID CIL-DLM-6196-10X-1.2
Unit 1,2 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Methylphosphonic acid (methyl-D3,98%)	[104801-16-3]	1000	µg/mL		