

PCB MIXTURE (34, 72, 100, 209) (unlabelled), 400 µg/mL in Isooctane

Art. ID CIL-EC-5605
Unit 1,2 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1262 II E2

Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2',3,5-Trichlorobipheny I (PCB 34)	[37680-68-5]	400	ng/mL		
	2,3',5,5'-Tetrachlorobi phenyl (PCB 72)	[41464-42-0]	400	ng/mL		
	2,2',4,4',6-Pentachloro biphenyl (PCB 100)	[39485-83-1]	400	ng/mL		
	Decachlorobiphenyl (PCB 209)	[2051-24-3]	400	ng/mL		