

EXTENDED PHENOLIC CALIBRATION SOLUTION CS3 ((unlabelled)/13C, 99%)

Art. ID CIL-ES-5510-CS3
Unit 0,5 mL
Deliverydetails De Minimis

Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
CS3	4-Methoxy-2,3,3',4',5-p entachlorobiphenyl		2,6	ng/mL		
CS3	4-Methoxy-2,2',3,4',5,5 '-hexachlorobiphenyl		2,59	ng/mL		
CS3	4-Methoxy-2,2',3,4',5,5 '-heptachlorobiphenyl		2,59	ng/mL		
CS3	5-Chloro-2-(2,4-dichlor ophenoxy)-anisole (meth yltriclosan)		26,2	ng/mL		
CS3	4-Methoxy-2,2',3,4'-tet rabromodiphenyl ether	[n/a]	2,57	ng/mL		
CS3	3-Methoxy-2,2',4,4'-tet rabromodiphenyl ether	[n/a]	2,57	ng/mL		
CS3	5-Methoxy-2,2',4,4'-tet rabromodiphenyl ether	[n/a]	2,57	ng/mL		
CS3	6-Methoxy-2,2',4,4'-tet rabromodiphenyl ether	[n/a]	2,57	ng/mL		
CS3	4'-Methoxy-2,2',4,5'-te trabromodiphenyl ether	[n/a]	2,57	ng/mL		
CS3	5'-Methoxy-2,2',4,4',5- pentabromodiphenyl ethe r	[n/a]	2,56	ng/mL		
CS3	6-Methoxy-2,2',4,4',5-P entabromodiphenyl ether	[497106-83-9]	2,56	ng/mL		
CS3	2,4-Dibromoanisole	[21702-84-1]	2,64	ng/mL		
CS3	2,4,5-Tribromoanisole	[n/a]	2,61	ng/mL		
CS3	2,4,6-Tribromoanisole	[607-99-8]	2,61	ng/mL		
CS3	Pentabromoanisole	[n/a]	2,57	ng/mL		
CS3	Pentachloroanisole	[1825-21-4]	26,3	ng/mL		

CS3	Dimethyl tetrabromobisp henol A	[n/a]	1,31	ng/mL
CS3	4-Methoxy-2,3,3',4',5-p entachlorobiphenyl (13C 12,99%)		104	ng/mL
CS3	4-Methoxy-2,2',3,4',5,5 '-hexachlorobiphenyl (1 3C12,99%)		104	ng/mL
CS3	4-Methoxy-2,2',3,4',5,5 '-6-heptachlorobiphenyl (13C12,99%)		103	ng/mL
CS3	5-Chloro-2-(2,4-dichlor ophenoxy)-anisole (meth yltriclosan) (ring-13C1 2,99%)		1047	ng/mL
CS3	6-Methoxy-2,2',4,4'-Tet rabromodiphenyl ether (ring-13C12,99%)	[n/a]	103	ng/mL
CS3	6'-Methoxy-2,2',4,4',6- pentabromodiphenyl ethe r (13C12,99%)		102	ng/mL
CS3	2,4-Dibromoanisole (rin g-13C6,99%)	[n/a]	105	ng/mL
CS3	2,4,5-Tribromoanisole (ring-13C6,99%)		104	ng/mL
CS3	2,4,6-Tribromoanisole (ring-13C6,99%)	[1097192-98-7]	104	ng/mL
CS3	Pentabromoanisole (13C6 ,99%)		103	ng/mL
CS3	Pentachloroanisole (13C 6,99%)	[1825-21-4]/Unl abeled	1052	ng/mL
CS3	Dimethyl Tetrabromobisp henol A (ring-13C12,99%)	[n/a]	53	ng/mL
CS3	1,2,3,4-Tetrachlorodibe nzo-p-dioxin (13C6,99%)	[114423-99-3]	25	ng/mL
CS3, IUPAC# 208	2,2',3,3',4,5,5',6,6'-N onachlorobiphenyl (13C1	[234432-92-9]	100	ng/mL

	2,99%) (IUPAC #208)			
CS3, IUPAC# 77	3,3',4,4'-Tetrabromodiphenyl ether (13C12,99%) (IUPAC #77)		75	ng/mL
CS3, IUPAC# 139	2,2',3,4,4',6-Hexabromodiphenyl ether (13C12,99%) (IUPAC #139)	[446254-96-2]/unlabelled	75	ng/mL