

JECS PHENOL/PARABEN SYRINGE STANDARD (D, 98%) in Methanol-OD

Art. ID	CIL-ES-5630
Unit	1,2 mL
Deliverydetails	ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Methyl paraben (methyl 4-hydroxybenzoate) (2,3,5,6-D4,98%)		250	ng/mL		
	n-Propyl paraben (n-propyl 4-hydroxybenzoate) (2,3,5,6-D4,98%)		250	ng/mL		
	p-n-Nonylphenol (ring-D4,OD,98%)		2500	ng/mL		
	Bisphenol S (ring-D8,98%)	[80-09-1]/unlabelled	250	ng/mL		
	Bisphenol F (D10,98%)	[1794786-93-8]	2500	ng/mL		
	Bisphenol A (3,3',5,5'-D4,97%)		2500	ng/mL		