

## **METHOD 8327 TARGET ANALYTE MIXTURE in Methanol (W/4 MOLAR EQUIV. NAOH)**

Art. ID CIL-ES-5642  
Unit 1,2 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

### Description

Gravimetric Concentration is based on the free acid, where applicable

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Perfluorobutyric acid (PFBA)	[375-22-4]	2000	ng/mL		
	Perfluoropentanoic acid (PFPeA), sodium salt	[2706-89-0]	2000	ng/mL		
	Perfluorohexanoic acid (PFHxA), sodium salt	[2923-26-4]	2000	ng/mL		
	Perfluoroheptanoic acid (PFHpA)	[375-85-9]	2000	ng/mL		
	Perfluoro-n-octanoic acid (PFOA)	[335-67-1]	2000	ng/mL		
	Perfluorononanoic acid (PFNA)	[375-95-1]	2000	ng/mL		
	Perfluorodecanoic acid (PFDA), sodium salt	[3830-45-3]	2000	ng/mL		
	Perfluoroundecanoic acid (PFUA), sodium salt	[60871-96-7]	2000	ng/mL		
	Perfluorododecanoic acid (PFDoA), sodium salt	[307-67-5]	2000	ng/mL		
	Perfluorotridecanoic acid (PFTrDA)	[72629-94-8]	2000	ng/mL		
	Perfluorotetradecanoic acid (PFTeDA)	[376-06-7]	2000	ng/mL		
	Perfluorobutanesulfonate (PFBS), potassium salt	[29420-49-3]	2000	ng/mL		
	Perfluoropentanesulfonate (PFPeS), sodium salt	[630402-22-1]	2000	ng/mL		
	Potassium perfluorohexanesulfonate (PFHxS) (mi	[3871-99-6]	2000	ng/mL		

x of isomers)				
Perfluoroheptanesulfonate (PFHpS), sodium salt	[21934-50-9]	2000		ng/mL
Perfluorooctanesulfonic acid (PFOS) (mix of isomers)	[1763-23-1]	2000		ng/mL
Perfluorononanesulfonate (PFNS), sodium salt	[98789-57-2]	2000		ng/mL
Perfluorodecanesulfonate (PFDS), sodium salt	[2806-15-7]	2000		ng/mL
1H,1H,2H,2H-Perfluorooxane sulfonate (4:2 FTS), sodium salt	[27619-93-8]	2000		ng/mL
1H,1H,2H,2H-Perfluorooxane sulfonate (6:2 FTS), sodium salt	[27619-94-9]	2000		ng/mL
1H,1H,2H,2H-Perfluorodecane sulfonate (8:2 FTS), sodium salt	[27619-96-1]	2000		ng/mL
Perfluorooctanesulfonamide (PFOSA)	[754-91-6]	2000		ng/mL
N-Methylperfluorooctane sulfonamidoacetic acid (N-MeFOSAA) (mix of isomers)	[2355-31-9]	2000		ng/mL
N-Ethylperfluorooctane sulfonamidoacetic acid (N-EtFOSAA) (mix of isomers)	[2991-50-6]	2000		ng/mL