

TCA CYCLE STANDARD 1 MIX (unlabelled)

Art. ID CIL-MSK-TCA1-US-1

Unit 1 vial

Description

Composition of the stable isotope-labeled TCA mixtures. Reconstituting a given vial in 1 mL of solvent (e.g., water) will produce an equimolar concentration of 100 µM.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Fumaric acid, disodium salt		100	µM		
	DL-2-HYDROXYGLUTARIC ID, DISODIUM SALT		100	µM		
	alpha-Ketoglutaric acid, disodium salt		100	µM		
	Malic acid, disodium salt		100	µM		
	SODIUM L-LACTATE		100	µM		
	Sodium pyruvate	[113-24-6]	100	µM		
	Succinic acid, disodium salt		100	µM		
	Trisodium citrate		100	µM		