

## **2,3,4',6-Tetrachlorobiphenyl 100 µg/mL IN ISOOCTANE (unlabelled), Certified Standard**

Art. ID	CIL-PCB-64-CS
Unit	1,2 mL
Deliverydetails	ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1262 II E2

### Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2,3,4',6-Tetrachlorobiphenyl (PCB 64)	[52663-58-8]	100	µg/mL		