

### **Custom mixture - 8 components 10 µg/ml each of in Methanol**

Art. ID LM24-CUS-10223  
Unit 2 x 10 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

#### Description

Haltbarkeit 18 Monate (original verschlossen)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	4-Heptylphenol	[1987-50-4]	10	µg/mL		
	4-n-Nonylphenol	[104-40-5]	10	µg/mL		
	4-n-Octylphenol	[1806-26-4]	10	µg/mL		
	4-n-Amylphenol	[14938-35-3]	10	µg/mL		
	4-tert-Amylphenol	[80-46-6]	10	µg/mL		
	4-tert-Butylphenol	[98-54-4]	10	µg/mL		
	4-tert-Octylphenol	[140-66-9]	10	µg/mL		
	2,4,6-Tri-tert-butylphe nol	[732-26-3]	10	µg/mL		