

Custom Solution - Phenole - Mix 4 100 µg/mL each of in Acetone

Art. ID LM24-CUS-10124-7X1
Unit 7 x 1 mL
Deliverydetails De Minimis UN1090 II E2

Description

Conservation: 3 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Pentachlorophenol acetate	[1441-02-7]	100	µg/mL		
Pentachlorophenylpyroglutamamat	L-Pyroglutamic Acid Pentachlorophenyl Ester	[28990-85-4]	100	µg/mL		
	Pentachlorophenol	[87-86-5]	100	µg/mL		
Pentachlorophenyldodecanoat	Pentachlorophenol dodecanoate	[3772-94-9]	100	µg/mL		