

**Custom Solution - 3 components 2000 mg/L each of in Acetonitrile / Dichloromethane tr.
(1)**

Art. ID LM24-CUS-71454
Unit 3 x 1 mL ampoule
Deliverydetails De Minimis UN1648 II E2

Description

Conservation: 12 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Benzidine-d8	[92890-63-6]	2000	mg/L		
	Anthracene D10	[1719-06-8]	2000	mg/L		
	2,4,5-Trichloroaniline	[636-30-6]	2000	mg/L		