

Reagecon Polychlorinated Biphenyl (PCB) Deuterated Internal Standard (2 Compound Mix) in Hexane

Art. ID REA-REPCB005-I
Unit 1 mL
Deliverydetails De Minimis UN1208 II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	3,3',4,4'-Tetrachlorobiphenyl-D6	[93952-23-9]	10	µg/mL		
	2,3,3',4,4',5-Hexachlorobiphenyl-2',6,6'-D3	[1219798-44-3]	10	µg/mL		