

COVID-19 Healgen Rapid Antigen Test Cassette Swab, 15 min Nasal Swab

Art. ID	REA-GCCOV-502a-H1
Unit	each
Deliverydetails	No Dangerous Good

Description

The Healgen Rapid COVID-19 Antigen Test Single Cassette, is an in vitro immunochromatographic assay for the qualitative detection of nucleocapsid protein antigen from SARS-CoV-2 in direct nasal swab from individuals who are suspected of COVID-19. The Healgen Rapid Covid-19 Antigen Test is intended to aid in the rapid diagnosis of SARS-CoV-2 infections (15 min Nasal Swab). This antigen test can be used to detect both wild and mutant SARS-CoV-2 variants of concern, including Omicron (B.1.1.529), Delta (B.1.617.2), UK (B.1.1.7), South Africa (B.1.351) and Brazil (P.1).