

Biphenyl 100 µg/mL in Acetonitrile

Art. ID LM24-N-15100-0010-100AN5
Unit 5 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1648 II E2

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|------------|---------------------|-------|--------|--------|
| | Biphenyl | [92-52-4] | 100 | µg/mL | | |