

Custom Solution - Chlorobenzenes - Mix 8 500 µg/ml each of in n-Hexane

Art. ID LM24-CUS-10225
Unit 5 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1208 II E2

Description

Conservation: 18 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	1,2,3-Trichlorobenzene	[87-61-6]	500	µg/mL		
	1,2,4-Trichlorobenzene	[120-82-1]	500	µg/mL		
	1,3,5-Trichlorobenzene	[108-70-3]	500	µg/mL		
	1,2,3,4-Tetrachlorobenzene	[634-66-2]	500	µg/mL		
	1,2,3,5-Tetrachlorobenzene	[634-90-2]	500	µg/mL		
	1,2,4,5-Tetrachlorobenzene	[95-94-3]	500	µg/mL		
	Pentachlorobenzene	[608-93-5]	500	µg/mL		
	Hexachlorobenzene	[118-74-1]	500	µg/mL		