

**Set of 2x10 mL Potassium Dichromate Solution for Absorbance Control 235 - 350 nm and 6x10 mL blank acc. 2.2.25 of Ph Eur**

Art. ID LM24-OP-4011669  
Unit 80 mL  
Deliverydetails No Dangerous Good /not restricted

Description

57 - 63 mg/l Potassium dichromate R in 0.005 M sulphuric acid and blank of 0.005 M sulphuric acid for the control of absorbance 235 - 350 nm /// This product conforms to the specifications given in 2.2.25 of the European Pharmacopoeia (old edition)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Potassium dichromate	[7778-50-9]	60,01 ± 0,15	mg/L		