

Set of 2x10 mL Potassium Dichromate Solution for Absorbance Control at 430 nm and 6x10 mL blank acc. 2.2.25 of Ph Eur

Art. ID LM24-OP-4011670
Unit 80 mL
Deliverydetails No Dangerous Good /not restricted

Description

570 - 630 mg/l Potassium dichromate R in 0.005 M sulphuric acid and blank of 0.005 M sulphuric acid for the control of absorbance at 430 nm /// This product conforms to the specifications given in 2.2.25 of the old European Pharmacopoeia.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Potassium dichromate	[7778-50-9]	600,4 ± 1,1	mg/L		