

Custom Solution - 8 component in HNO3 5%

Art. ID	LM24-CUS-71477
Unit	1 x 125 mL bottle
Deliverydetails	ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Description

Conservation 24 months (originally closed bottle), 12 months after opening.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Antimony (Sb)	[7440-36-0]	500	mg/L		
	Boron (B)	[7440-42-8]	500	mg/L		
	Arsenic (As)	[7440-38-2]	500	mg/L		
	Selenium (Se)	[7782-49-2]	500	mg/L		
	Thallium (Tl)	[7440-28-0]	500	mg/L		
	Manganese (Mn)	[7439-96-5]	1000	mg/L		
	Tin (Sn)	[7440-31-5]	1000	mg/L		
	Phosphorus (P)	[7723-14-0]	2000	mg/L		