

EPA Method 525 Internal Standards PAH Mixture 2000 µg/mL in Dichloromethane

Art. ID DRE-A50000277DI
Unit 1 mL
Deliverydetails De Minimis UN1593 III E1

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|--------------|---------------------|------|--------|--------|
| | Acenaphthene D10 | [15067-26-2] | | | | |
| | Chrysene D12 | [1719-03-5] | | | | |
| | Perylene D12 | [1520-96-3] | | | | |
| | Phenanthrene D10 | [1517-22-2] | | | | |