

**EPA Method 525 Internal Standards PAH Mixture 2000 µg/mL in Dichloromethane**

Art. ID DRE-A50000277DI  
Unit 1 mL  
Deliverydetails De Minimis UN1593 III E1

| Text/Information | Analyte/Parameter | CAS number   | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|--------------|---------------------|------|--------|--------|
|                  | Acenaphthene D10  | [15067-26-2] |                     |      |        |        |
|                  | Chrysene D12      | [1719-03-5]  |                     |      |        |        |
|                  | Perylene D12      | [1520-96-3]  |                     |      |        |        |
|                  | Phenanthrene D10  | [1517-22-2]  |                     |      |        |        |