

Copper base, Type: Brass, RM, Mushroom 60 mm x 5 mm, batch 2371-223

Art. ID CTIF-UZ51-2371-223
Unit each
Deliverydetails No Dangerous Good /not restricted

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Arsenic (As) | [7440-38-2] | 0,085 | % | | |
| | Copper (Cu) | [7440-50-8] | 83,82 | % | | |
| | Iron (Fe) | [7439-89-6] | 0,10 | % | | |
| | Manganese (Mn) | [7439-96-5] | <0,001 | % | | |
| | Nickel (Ni) | [7440-02-0] | 0,155 | % | | |
| | Phosphorus (P) | [7723-14-0] | 0,01 | % | | |
| | Lead (Pb) | [7439-92-1] | 0,19 | % | | |
| | Silicon (Si) | [7440-21-3] | 0,01 | % | | |
| | Tin (Sn) | [7440-31-5] | 1,52 | % | | |
| | Zinc (Zn) | [7440-66-6] | 14,1 | % | | |