

Custom Solution - 4 components 200 µg/mL each of in Methanol

Art. ID LM24-CUS-71513
Unit 3 x 1 mL ampoule
Deliverydetails De Minimis UN1230 (6.1) II E2

Description

Conservation 12 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Epichlorhydrin	[106-89-8]	200	µg/mL		
	1,3-Dichloropropan-2-ol	[96-23-1]	200	µg/mL		
	1,2,3-Trichloropropane	[96-18-4]	200	µg/mL		
	2,3-Dichloro-1-propanol	[616-23-9]	200	µg/mL		