

Iron Base, Steel - Stainless, SUS, Set of KTC-5-31 - KTC-5-38 only, Batch B, see KTC-5-31-38-SET, Disc 35 mm x 20 mm

Art. ID KTC-5-36
Unit disc
Deliverydetails No Dangerous Good /not restricted

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Soluble Al | | 0,001 | % | | |
| | Insoluble Al | | 0,008 | % | | |
| | Arsenic (As) | [7440-38-2] | <0,001 | % | | |
| | Carbon (C) | [7440-44-0] | 0,15 | % | | |
| | Cobalt (Co) | [7440-48-4] | 0,003 | % | | |
| | Chromium (Cr) | [7440-47-3] | 22,31 | % | | |
| | Copper (Cu) | [7440-50-8] | <0,01 | % | | |
| | Manganese (Mn) | [7439-96-5] | 0,43 | % | | |
| | Molybdenum (Mo) | [7439-98-7] | 0,043 | % | | |
| | Niobium (Nb) | [7440-03-1] | 0,001 | % | | |
| | Nickel (Ni) | [7440-02-0] | 0,11 | % | | |
| | Phosphorus (P) | [7723-14-0] | 0,014 | % | | |
| | Sulfur (S) | [7704-34-9] | 0,009 | % | | |
| | Silicon (Si) | [7440-21-3] | 0,14 | % | | |
| | Titanium (Ti) | [7440-32-6] | 0,005 | % | | |