

**see VASKUT-T3/2-40MM - Iron Base, Steel - Low Alloy, CRM, Transformer, Dynamo, Disc
32-35 mm x 39 mm, low stock - please ask for availability**

Art. ID VASKUT-T3/2
Unit disc
Deliverydetails No Dangerous Good /not restricted

Description

Also available in 30-35 mm x 19 mm - please ask for further information.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Carbon (C)	[7440-44-0]	0,09	%		
	Chromium (Cr)	[7440-47-3]	0,4	%		
	Copper (Cu)	[7440-50-8]	0,1	%		
	Manganese (Mn)	[7439-96-5]	0,6	%		
	Nickel (Ni)	[7440-02-0]	0,11	%		
	Phosphorus (P)	[7723-14-0]	0,058	%		
	Sulfur (S)	[7704-34-9]	0,033	%		
	Silicon (Si)	[7440-21-3]	0,66	%		
informational value	Titanium (Ti)	[7440-32-6]	<0,01	%		