

Quality Control Standard of Mineral Oils (ISO 9377-2-Mod) - 2 components 10000 mg/L each of in n-Hexane

Art. ID LM24-OP-4011724
Unit 10 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1208 II E2

Description

Validity: 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Mineral Oil (CAS: [n/a]))	[n/a]	10000	mg/L		
	Diesel Oil (CAS: [n/a])	[n/a]	10000	mg/L		