

**Quality Control Standard of Mineral Oils (ISO 9377-2-Mod) - 2 components 10000 mg/L each of in n-Hexane**

Art. ID LM24-OP-4011724  
Unit 10 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1208 II E2

Description

Validity: 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Mineral Oil [CAS:74869-22-0]	[74869-22-0]	10000	mg/L		
	Diesel Oil [CAS:68334-30-5]	[68334-30-5]	10000	mg/L		