

n-Alkane Standard Solution 16 components C10-C40 (all even) (ISO 9377-2-Mod) 100 mg/L each of in n-Hexane / Petroleum ether (1/1)

Art. ID LM24-OP-4011727
Unit 1 mL
Deliverydetails De Minimis UN1993 II E2

Description

Validity: 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	n-Decane	[124-18-5]	100	mg/L		
	n-Dodecane	[112-40-3]	100	mg/L		
	n-Tetradecane	[629-59-4]	100	mg/L		
	n-Hexadecane	[544-76-3]	100	mg/L		
	n-Octadecane	[593-45-3]	100	mg/L		
	n-Eicosane	[112-95-8]	100	mg/L		
	n-Docosane	[629-97-0]	100	mg/L		
	n-Tetracosane	[646-31-1]	100	mg/L		
	n-Hexacosane	[630-01-3]	100	mg/L		
	n-Octacosane	[630-02-4]	100	mg/L		
	n-Triacontane	[638-68-6]	100	mg/L		
	n-Dotriacontane	[544-85-4]	100	mg/L		
	n-Tetratriacontane	[14167-59-0]	100	mg/L		
	n-Hexatriacontane	[630-06-8]	100	mg/L		
	n-Octatriacontane	[7194-85-6]	100	mg/L		
	n-Tetracontane	[4181-95-7]	100	mg/L		