

## **6 components 100 µg/mL each of in Methanol**

Art. ID                      LM24-OP-4011824  
Unit                         1 mL  
Deliverydetails            De Minimis UN1230 (6.1) II E2

### Description

Validity: 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Diuron	[330-54-1]	100	µg/mL		
	Isoproturon	[34123-59-6]	100	µg/mL		
	Simazine	[122-34-9]	100	µg/mL		
	Terbutryn	[886-50-0]	100	µg/mL		
	Irgarol 1051	[28159-98-0]	100	µg/mL		
	Atrazine	[1912-24-9]	100	µg/mL		