

8 components 100 µg/mL each of in Acetonitrile

Art. ID LM24-OP-4011840
Unit 1 mL
Deliverydetails De Minimis UN1648 II E2

Description

Validity: 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Diuron	[330-54-1]	100	µg/mL		
	Isoproturon	[34123-59-6]	100	µg/mL		
	Prometryn	[7287-19-6]	100	µg/mL		
	Propazine	[139-40-2]	100	µg/mL		
	Simazine	[122-34-9]	100	µg/mL		
	Terbutylazine	[5915-41-3]	100	µg/mL		
	Terbutryn	[886-50-0]	100	µg/mL		
	Atrazine	[1912-24-9]	100	µg/mL		