

8 components 100 µg/mL each of in Methanol

Art. ID LM24-OP-4012113
Unit 1 mL
Deliverydetails De Minimis UN1230 (6.1) II E2

Description

Validity: 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Chloroform	[67-66-3]	100	µg/mL		
	1,2-Dichloroethane	[107-06-2]	100	µg/mL		
	Ethylbenzene	[100-41-4]	100	µg/mL		
	Toluene	[108-88-3]	100	µg/mL		
	m-Xylene	[108-38-3]	100	µg/mL		
	o-Xylene	[95-47-6]	100	µg/mL		
	p-Xylene	[106-42-3]	100	µg/mL		
	Benzene	[71-43-2]	100	µg/mL		