

18 components 200 µg/mL each of in Methanol

Art. ID LM24-OP-4012131
Unit 1 mL
Deliverydetails De Minimis UN1230 (6.1) II E2

Description

Validity: 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Chloroform	[67-66-3]	200	µg/mL		
	cis-1,2-Dichloroethene	[156-59-2]	200	µg/mL		
	Chlorobenzene	[108-90-7]	200	µg/mL		
	Dibromochloromethane	[124-48-1]	200	µg/mL		
	1,2-Dichloroethane	[107-06-2]	200	µg/mL		
	Ethylbenzene	[100-41-4]	200	µg/mL		
	Styrene	[100-42-5]	200	µg/mL		
	Tetrachloroethene	[127-18-4]	200	µg/mL		
	trans-1,2-Dichloroethene	[156-60-5]	200	µg/mL		
	Bromoform (Tribromomethane)	[75-25-2]	200	µg/mL		
	Trichloroethene	[79-01-6]	200	µg/mL		
	Toluene	[108-88-3]	200	µg/mL		
	Vinylchloride	[75-01-4]	200	µg/mL		
	m-Xylene	[108-38-3]	200	µg/mL		
	p-Xylene	[106-42-3]	200	µg/mL		
	1,1-Dichloroethene	[75-35-4]	200	µg/mL		
	Bromodichloromethane	[75-27-4]	200	µg/mL		
	Benzene	[71-43-2]	200	µg/mL		