

Iron Base, Steel - Stainless, SUS, 304L, Disc 40 mm x 25 mm

Art. ID TEC-PV-111/1
Unit disc
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Carbon (C)	[7440-44-0]	0,0226	%		
	Chromium (Cr)	[7440-47-3]	18,49	%		
	Copper (Cu)	[7440-50-8]	0,105	%		
	Manganese (Mn)	[7439-96-5]	1,538	%		
	Molybdenum (Mo)	[7439-98-7]	0,173	%		
	Nickel (Ni)	[7440-02-0]	8,57	%		
	Phosphorus (P)	[7723-14-0]	0,019	%		
	Sulfur (S)	[7704-34-9]	0,026	%		
	Silicon (Si)	[7440-21-3]	0,485	%		