

### **Custom Solution - 6 components 2000 µg/mL each of in Methanaol**

Art. ID LM24-CUS-71517  
Unit 4 x 1 mL ampoule  
Deliverydetails De Minimis UN1230 (6.1) II E2

#### Description

Conservation 12 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Chloroform (Trichlorome thane)	[67-66-3]	2000	µg/mL		
	Bromoform (Tribromometh ane)	[75-25-2]	2000	µg/mL		
	Bromodichloromethane	[75-27-4]	2000	µg/mL		
	Dibromochloromethane	[124-48-1]	2000	µg/mL		
	Tetrachloroethene	[127-18-4]	2000	µg/mL		
	Trichloroethene	[79-01-6]	2000	µg/mL		