

Cast Iron and Pig Iron, CRM, truncated cone h 38-42mm, Ø of top base- 36-40 mm, of bottom base- 38-42mm

Art. ID VS-ChG10/5
Unit each
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silicon (Si)	[7440-21-3]	2,79	%		
	Manganese (Mn)	[7439-96-5]	0,89	%		
	Sulfur (S)	[7704-34-9]	0,11	%		
	Phosphorus (P)	[7723-14-0]	0,095	%		
base	Iron (Fe)	[7439-89-6]				