

**C20 Sphingomyelin (d18:1/20:0)**

Art. ID LA-56-4020-45  
Unit 500 µg

**Description**

Lipid number: d18:1/20:0 /// Documentation: Certificate of Analysis /// Purity: >98%

| Text/Information | Analyte/Parameter              | CAS number    | Concentration/Value | Unit | Method | Source |
|------------------|--------------------------------|---------------|---------------------|------|--------|--------|
|                  | C20 Sphingomyelin (d18:1/20:0) | [121999-68-6] |                     |      |        |        |