

Custom solution - Sorbic acid [110-44-1] 10000 µg/mL in Methanol

Art. ID LM24-CUS-71537
Unit 3 x 1 mL
Deliverydetails De Minimis UN1230 (6.1) II E2

Description

Conservation: 12 months (originally closed bottle).

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|------------|---------------------|-------|--------|--------|
| | Sorbic acid | [110-44-1] | 10000 | µg/mL | | |