

Custom Solution - 18 components in HNO3 2%

Art. ID LM24-CUS-71559
Unit 1 x 500 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

Description

Conservation 8 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	40	mg/L		
	Arsenic (As)	[7440-38-2]	0,2	mg/L		
	Boron (B)	[7440-42-8]	0,8	mg/L		
	Calcium (Ca)	[7440-70-2]	100	mg/L		
	Cobalt (Co)	[7440-48-4]	0,2	mg/L		
	Iron (Fe)	[7439-89-6]	40	mg/L		
	Manganese (Mn)	[7439-96-5]	2	mg/L		
	Molybdenum (Mo)	[7439-98-7]	0,2	mg/L		
	Sodium (Na)	[7440-23-5]	40	mg/L		
	Sulfur (S)	[7704-34-9]	40	mg/L		
	Antimony (Sb)	[7440-36-0]	0,2	mg/L		
	Selenium (Se)	[7782-49-2]	0,2	mg/L		
	Tin (Sn)	[7440-31-5]	0,2	mg/L		
	Tellurium (Te)	[13494-80-9]	0,1	mg/L		
	Thallium (Tl)	[7440-28-0]	0,2	mg/L		
	Vanadium (V)	[7440-62-2]	0,2	mg/L		
	Tungsten (W)	[7440-33-7]	0,1	mg/L		
	Beryllium (Be)	[7440-41-7]	0,2	mg/L		