

Copper Alloy, Leaded Red Brass, CRM, CDA 836 / C83600, chips

Art. ID IARM-CU836-18-C
Unit 100 g (chips)
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Antimony (Sb)	[7440-36-0]	0,103	%		
	Arsenic (As)	[7440-38-2]	0,0098	%		
	Bismuth (Bi)	[7440-69-9]	0,035	%		
	Cadmium (Cd)	[7440-43-9]	0,0014	%		
	Cobalt (Co)	[7440-48-4]	0,0026	%		
	Copper (Cu)	[7440-50-8]	84,5	%		
	Iron (Fe)	[7439-89-6]	0,049	%		
	Lead (Pb)	[7439-92-1]	4,9	%		
	Nickel (Ni)	[7440-02-0]	0,423	%		
	Niobium (Nb)	[7440-03-1]	0,002	%		
	Phosphorus (P)	[7723-14-0]	0,028	%		
	Selenium (Se)	[7782-49-2]	0,0016	%		
	Silver (Ag)	[7440-22-4]	0,0292	%		
	Sulfur (S)	[7704-34-9]	0,041	%		
	Tin (Sn)	[7440-31-5]	5,02	%		
	Zinc (Zn)	[7440-66-6]	5,09	%		