

### **Custom Solution - 8 components 100 mg/L each in Methanol**

Art. ID LM24-CUS-71567  
Unit 5 x 1 mL ampoule  
Deliverydetails No Dangerous Good UN1230 (6.1) II E2

#### Description

Conservation: 12 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Bisphenol A	[80-05-7]	100	mg/L		
	Bisphenol B	[77-40-7]	100	mg/L		
	Bis(4-hydroxyphenyl)met hane	[620-92-8]	100	mg/L		
	alpha,alpha'-Bis(4-hydr oxyphenyl)-1,4-diisopro pylbenzene	[2167-51-3]	100	mg/L		
	4,4'-Sulfonyldiphenol	[80-09-1]	100	mg/L		
	1,1-Bis(4-hydroxyphenyl )cyclohexane	[843-55-0]	100	mg/L		
	2,2-Bis(4-hydroxyphenyl )hexafluoropropane	[1478-61-1]	100	mg/L		
	4,4'-(alpha-Methylbenzy lidene)bisphenol	[1571-75-1]	100	mg/L		