

Diesel oil (type A) (ISO 9377-2, ISO 16703, EN 14039)

Art. ID	BAM-K008a
Unit	2 mL
Deliverydetails	No Dangerous Good /not restricted

Description

Diesel oil, certified mass fraction of the boiling range C10 - C40 to be used as calibration standard (type A) for the determination of TPH by GC/FID. Certified value

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Mass fraction of the boiling range C10 - C40		0,942 ± 0,013	g/g		