

Lubricating oil (type B) (ISO 9377-2, ISO 16703, EN 14039)

Art. ID BAM-K009a
Unit 2 mL
Deliverydetails No Dangerous Good /not restricted

Description

Lubricating oil, certified mass fraction of the boiling range C10 - C40 to be used as calibration standard (type B) for the determination of TPH by GC/FID and recovery standard for the gravimetric TPH determination. Certified value

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Mass fraction of the boiling range C10 - C40		0,982 ± 0,011	g/g		