

PFT - Mix 10 2 µg/ml each

Art. ID LM24-CUS-10231
Unit 1 x 10 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

Description

Haltbarkeit: 18 Monate (original verschlossen)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Perfluorooctane sulfona mide PFOSA	[754-91-6]	2	µg/mL		
	N-Methylperfluorooctane sulfonamide	[31506-32-8]	2	µg/mL		
	Sulfuramid	[4151-50-2]	2	µg/mL		
Perfluor-3,7-dimethyl- tansäure	2,2,3,4,4,5,5,6,6,7,8,8, ,8-Tridecafluoro-3,7-bi s(trifluoromethyl)-octa noic Acid	[172155-07-6]	2	µg/mL		
	7H-Perfluoroheptanoic a cid	[1546-95-8]	2	µg/mL		
	4,4,5,5,6,6,7,7,8,8,9,9, ,10,10,11,11,11-Heptade cafluoroundecanoic acid	[34598-33-9]	2	µg/mL		
	1H,1H,2H,2H-Perfluoro-n -octane sulfonic acid	[27619-97-2]	2	µg/mL		
	Perfluoroheptanesulfoni c Acid	[375-92-8]	2	µg/mL		
	1H,1H,2H,2H-Perfluorode canesulfonic acid	[39108-34-4]	2	µg/mL		
	1H,1H,2H,2H-Perfluorodo decane sulfonic acid	[120226-60-0]	2	µg/mL		