

Cold Filter Plugging Point (CFPP) Standard acc to ASTM D6371 and EN 116 (summer diesel)

Art. ID LM24-OP-4012526
Unit 250 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1202 III E2

Description

Validity: 24 months

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-----------------------------------|------------|---------------------|------|------------------------|--------|
| | Cold filter plugging point (CFPP) | | (-3,0 ± 0,1) | °C | ASTM D6371 / EN 116 | |